

# A Rapid and Scalable 96-Well Platform to Screen CARs and TCRs

Alexander Boucher<sup>1</sup>, Peter Moua<sup>1</sup>, Elise Kim<sup>1</sup>, Chris Driscoll<sup>1</sup>, Shen Lu<sup>1</sup>, Joshua Ferrell<sup>2</sup>, Andrew Basinski<sup>3</sup>, Rukmini Ladi<sup>4</sup>, Stacie Seidel<sup>1</sup>, Cherylene A. Plewa<sup>1</sup> and Ana Maria Chadbourne<sup>1</sup>  
<sup>1</sup>ElevateBio, 225 Wyman St, Waltham MA; <sup>2</sup>Current Affiliation: Affini-T Therapeutics, Watertown, MA; <sup>3</sup>Current Affiliation: OneCyte Biotechnologies, Cambridge, MA; <sup>4</sup>Current Affiliation: Sartorius Stedim Biotech


## Abstract

- Identification of a binder sequence that will give rise to a chimeric antigen receptor (CAR) or T cell receptor (TCR) construct is a necessary first step in the development of engineered cell therapies.
- Screening platforms for CAR and TCR construct nomination in primary T cells are seldom scalable, utilize non-MOI normalized lentiviral vector preparations and rely on limited phenotypical or functional assessments.
- Because of these limitations, the performance of top early candidates may not translate to larger-scale preparations, thus requiring time consuming repeat experiments to achieve reproducible data at-scale for *in vivo* work.
- Here we report a **rapid, mid-throughput, comprehensive, scalable CAR/TCR screening platform for the selection of candidate CARs and TCRs in primary T cells.**
- The platform relies on a **scalable mid-throughput lentiviral vector production platform**, and a **96-well scalable T cell expansion that produces enough cells for characterization and functional assays** including vector copy number (VCN), surface binder expression and function.
- The platform was utilized to **screen candidate CARs by surface expression and VCN, tonic signaling, functional activation by engineered cell lines as well as cell lines expressing endogenous high, medium and low levels of target.**
- Within **three rounds of screening**, we were able to nominate **top candidates ready for *in vivo* testing.**
- Our **proprietary CAR/TCR screening platform thus allows acceleration of nomination of binders for faster development of CAR and TCR therapies.**


## Summary of Screening Workflow

### Criteria for Candidate Selection:


1. CAR Expression at equal or higher than 30%
2. IFN $\gamma$  production and cytotoxic function equal or higher than control, lack of tonic signaling
3. Diversity in B-Cell lineage for top clones
4. Sensitivity to target level expression in cell lines




## Lentiviral Platform and Target Cell Lines




## Screens 2&3: Cell Lines Expressing Varying Levels of Target Enable Top CAR Sequence Selection




## Screen 1: CAR Expression, IFN $\gamma$ and Cytotoxicity Data for 73 CAR Sequences



## Screen 4: *in vitro* Validation of Top Three Candidates in Large Scale 100M G-Rex Shows Scalability, and Superior CAR Activity Compared to Control




## Small to Large Scale Activity Comparison of Top Three CAR Candidates



## Platform Validation for NYESO TCR

Transduction and Expansion Scalability was demonstrated with NY-ESO TCR as Benchmark



## Conclusion

- Our 96-well screening platform was validated for CARs & TCRs, and enabled **selection of candidate CARs within three rounds of screening.**
- The throughput and cell yield were sufficient to perform an extensive array of assays: **surface expression and VCN, tonic signaling, functional activation by engineered cell lines as well as endogenously expressing cell lines.**
- Inclusion of cell lines **expressing endogenous high, medium and low levels of target enabled higher stringency to select top candidates.**
- The **platform data is reproducible, and the platform is scalable.**
- Our **proprietary CAR/TCR platform thus allows acceleration of nomination of binders for faster development of CAR and TCR therapies.**